MITC (METHYL ISOTHIOCYANATE)



1. PERFORMANCE

1) Measuring range 0.66-22 ppm Number of pump strokes $1 (100 \text{m} \ell)$ $1/2 (50 \text{m} \ell)$ 1/2 minutes/1 pump stroke

3) Detectable limit : 0.1 ppm

4) Shelf life : 1 year (Necessary to store in refrigerated conditions ; $0 \sim 10^{\circ}$ C)

5) Operating temperature : $0 \sim 40 \,^{\circ}\text{C}$

6) Temperature compensation : Necessary (See "TEMPERATURE CORRECTION TABLE")
7) Reading : Direct reading from the scale calibrated by 1 pump stroke

8) Colour change : Pink → Yellowish orange

2. RELATIVE STANDARD DEVIATION

RSD-low: 10% RSD-mid.: 5% RSD-high: 5%

3. CHEMICAL REACTION

By decomposing with an oxidizer, Sulphur dioxide is produced and PH indicator is discoloured.

 $CH_3NCS + V_2O_5 + H_2SO_4 \rightarrow SO_2$

SO₂ + pH indicator (Pink) → Yellowish orange

4. CALIBRATION OF THE TUBE

DIFFUSION TUBE METHOD

5. INTERFERENCE AND CROSS SENSITIVITY

Substance	Interference	%	Coexistence
Carbon dioxide		35	Higher readings are given.

(NOTE)

The scale is calibrated based on the temperature of $20\,\mathrm{C}$. Reading obtained in other temperature circumstances should be corrected with the following temperature correction coefficient table.

TEMPERATURE CORRECTION COEFFICIENT TABLE(AT 20°C)

Pump stroke	Temperature (℃)	0	5	10	15	20	25	30	35	40
1	Coefficient	1.30	1.23	1.15	1.08	1.00	0.95	0.90	0.85	0.80
1/2	Coemcient	1.70	1.48	1.30	1.11	1.00	0.97	0.95	0.92	0.90

Actual concentration = Reading value × Coefficient for temperature correction